DENTAL CAD/CAM

GN-I *NOW Debut!*

"Dental CAD/CAM GN-I" is a computer assisted next-generation dental laboratory system that will bring a radical revolution to the conventional method of prosthetic preparation. Based on the fundamental study of Next-generation Oral Device Engineering System Project (presented by Nikon Corporation, Hitachi Seikoh Co., Ltd. and GC Corporation). During the period 1993 to 1997, Nikon developed the measuring device and layout designing software whilst GC developed the processing machine, prosthetic material, special-purpose bite and model materials. This system performs measurement, design and processing using the indirect method after taking impressions and study models of dentition with missing teeth.

Next-generation Oral Device Engineering System Project is a research project entrusted to Technology Research Association of Medical and Welfare Apparatus by New Energy and Industrial Technology Development Organization (NEDO), which was financed by Agency of Industrial Science and Technology, Ministry of International Trade and Industry in Japan.

#11, #21 restored with high compression composite resin block

#44, #45 restored with ceramic

#36 restored with titanium



Dental CAD/CAM system DENTAL CAD/CAM GN-I

System containing:
GN-I Measuring machine, GN-I Software,
GN-I Milling machine

The system does not contain a personal computer.
The system is produced on an indent basis with a delivery time of approximately twenty six weeks.

Precise measurement of undercut area by adoption of five axis measurement

"GN-I measuring machine" was designed for use in the laboratory using a laser sensor to determine the inner and outer contour of the restoration. By adopting a five axis measuring- model rotation (θ) and 180 degree rotation of the laser sensor (α) as well as length, width and height (x, y, z), the shape of undercut area can be replicated precisely. Measuring pitch is selective from 0.1, 0.3 and 0.5mm. Ordinal measurement is performed in 0.3mm pitch, and the minimum pitch is selected automatically for a margin in which accurate information is required. You can obtain precise data in a short time.


Five axis measurement- model rotation (θ) and 180 degree rotation of the laser sensor (α) as well as length, width and height (x, y, z)

The margins are set automatically. These can be easily adjusted using the computer mouse if any marginal areas are unclear.


Productive software function allowing free design

The restorative prosthesis is designed by enlargement, reduction and adjustment of one of the selected computer stored reference designs. The prosthetic design function is undertaken easily using a conventional mouse. If you input data you can adapt the software to store and create your own database. The software also makes it possible to copy the exact size shape and contour of the same tooth in the opposing arch to obtain a duplicate restoration.




The restoration can be easily viewed from any position by clicking on the appropriate design view button Example occlusal surface, distal, mesial, lingual and buccal.

Selecting the shading mode can highlight a curve with mesh, or make a model translucent as well as shadowed real image.




Prosthesis design and adjustment of shape and size is very similar to conventional laboratory waxing up

Design activities such as additional contouring of the proximal surface and adjustment of fissures are possible easily by selecting the size of the reference point (circle) and appointing the range. You can check the thickness simply because you can identify the sectional profile of a designed material.

Easy adjustment of contact

Occlusal surface contact and proximal tooth contact are established by overlapping measuring data obtained from special-purpose CAD silicone and proximal tooth measuring data respectively. The contact area on the abutment tooth is highlighted by masking the non contact areas with a different color.


Cement space to provide precision adaptability

Sufficient space for the luting can be programmed into computer design before commencing the milling procedure. Set the starting position and cement thickness required in proportion to the cement space and margin line, for high-precision prosthetic material bonding.

Automatic tool replacing function with superior working efficiency



Each milling tool is selected and changed automatically based on design. Even during production, you can independently design and preparing further restorations to be produced once the current job has been completed.


There are special GN-I milling machine burs for ceramic or composite restorations and five different milling burs for titanium preparation. These are contained in two separate holders with a duplicate set of burs holder.

Three different restorative materials to suit your particular design and functional requirements

Three different restorative materials are supplied for "GN-I". High Compression prefabricated composite resin and ceramic blocks providing superior flexural strength and wear and abrasion resistance. The third material is lightweight, strong, flexible, easy to polish titanium.


GN-I Measuring machine

(Measuring unit) outer diameter: 400×552×500mm, Weight: 42kg, Power: AC100V 50/60Hz, Object to measure: sectional and removable plaster model, bite impression of occlusal surface of antagonistic tooth, Measuring method: triangle measurement with laser light, Measuring accuracy: ±20μm, Measuring pitch: 0.1, 0.3, 0.5mm, Measuring time: three minutes per one tooth (in the case of measuring pitch:0.3mm), Measuring area: 90(x)×60(y)×25(z)mm (Power supply unit) : Outer diameter: 400×550×252mm, Weight: 38kg

GN-I Milling machine

Outer diameter: 610×610×1600mm, Weight: 200kg, Power: AC100V 50/60Hz, Material to process: GN-I high compression composite resin, GN-I ceramic block, GN-I titanium block, Axis transfer: 40(x)×30(y)×30(z)mm, 360 degree (A), Spindle rotation speed: Max 50,000rpm, Automatic tool replacement, Tool length compensation, Tool life management function,

GN-I Software

Prosthetic material can be designed for : inlays, crowns and copings.

Main functions: reference model registration, margin line extraction, proximal tooth contact setting, antagonistic tooth contact setting, cement space setting, shape adjustment, preparation of data for processing, data save.

PC configuration required: CPU: Pentium II, 266MHz, Hard disk capacity: 4GB, memory: 128MB, 3D graphic accelerator compliant to Open GL, Video memory: 8MB or more, Monitor: 1280×1024 resolution, LAN port, RS232-C port, CD-ROM drive, FD drive, OS : Windows NT4.0,

*Windows and Pentium II are the register trademarks of Microsoft corporation and Intel corporation.

*Specification and appearance of the product may be modified without notice.

GC CORPORATION

76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan Phone:03-3558-5181 Fax:09-3966-1470
<http://www.gcdental.co.jp/>

Prior to use, carefully read the instructions for use included in the package.

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Dear Dr. Duret,

Thank you very much for your recent letter . It has been already three month or so since we met you in Yokohama. I hoped you had a nice new year's day.

I will enclose some slides made for computer presentation. The file name is "CADCAMGN-1.ppt". It is recorded on CD enclosed.

If you have any questions on these slides, please feel free to ask me or Mr. Hirota, GC corporation.

Sincerely yours,

Makoto Ogino
System Integration Department
System Development Headquarters
NIKON CORPORATION
Tel:81-3-3773-5632
Fax: 81-3-3773-6357
ogino@nikongw.nikon.co.jp